



---

# Immersion OIL

---

**NEW  
PACKAGING**



*Comfort of use*

*Transparency*

*Cleanliness*

*User safety*

# What makes RAL immersion oil **UNIQUE ?**



## COMFORT OF USE

The profiled vial DropStick® has no cap and comes with a built-in dropper.

Thanks to the DropStand® rack you can always have the pair within easy reach.



## TRANSPARENCY

Through its high quality, RAL Immersion Oil helps optimize the reading of your slides.

RAL Immersion Oil has an unequalled level of transparency: it has the same refraction index as glass.



## CLEANLINESS

The viscosity of the oil has been designed to pour one drop at a time. The concept DropStick® - DropStand® allows an accurate use and no leak on your workbench.

**15 ml = up to 650 slides, approximately 650 drops**



## USER SAFETY

RAL Immersion Oil contains no toxic substances.

## RAL IMMERSION PACKS AVAILABLE



DropStand® with DropStick (1 x 15 ml + 1 DropStand) Ref.: 340060-0000

Pack DropStand® Immersion 30 (2 x 15 ml + 1 DropStand) Ref.: 340032-0000

Pack DropStand® Immersion 100 (7 x 15 ml + 1 DropStand) Ref.: 340103-0000

Pack Immersion 150 (6 x 25 ml + 2 built-in dropper) Ref.: 340150-0000

Pack Immersion 300 (6 x 50 ml + 2 built-in dropper) Ref.: 340300-0000

